


# Probe4Light

## Spectral Characterization of Light Sources



The lighting industry undergoes a major change with the introduction of the LED technology. The latter offers many advantages — ultra compact, low power consumption, optimize visual comfort — but involves to manage carefully its optical, electrical and thermal parameters.

Probe4Light delivers lab class spectra and colors measurements of LEDs. The spectrometer sensor — not RGB — enables to control accurately various types of LEDs like high power LEDs, monochromatic LEDs and polychromatic chips.

As a flexible hand held, standalone unit, Probe4Light also features the strength of a software interface to display all the device parameters and the color results at a glance.

The solution is suited for diagnostic, incoming inspection and quality control in laboratory as well as industrial and in-situ environments.

Although Probe4Light is initially designed for the new energy saving, LED-based products, it also works with other light source technologies.

Contents	
Measurement Configurations	Integrating sphere Diffuser (opal material) SMA-905 adapter
Sensor	Linear 1024 pxl CMOS sensor Diffraction grating True spectral, 380 — 780 nm range
Measurement Time	1 ms to 10 sec.
Interface	USB 2.0, self powered
Options	Other measurement heads (contact us)
Applications	Light source benchmark analysis LED characterization Lighting components incoming inspection

# Probe4Light


## Spectral Characterization of Light Sources

### Device description

- Stiff aluminium body with embedded technological blocks : spectrophotometer, electronic card and integrating sphere.
- 30 mm integrating sphere and 20 mm diffuser measurement heads. Other dimensions available on request. Inner sphere surface coated with barium sulphate.
- Lab class, high resolution spectrophotometer — 128 to 1024 pixels — optimized for spectral and colorimetric analysis. Match CIE Publ. No.15.2, Colorimetry, Second Edition (1986).
- Integration of time of the sensor is adjustable , manually or automatically — to fit better the measured light intensity. Also compatible with measurements of flashed or dimmed lamps.



- USB interface for driving the spectrophotometer and for reading the spectral data. No need for external power supply.
- PhotonProbe Software under Windows graphical interface. The user controls all the spectrophotometer parameters and manages to display the results in various numerical and diagram forms for a first analysis.
  - Spectrum 380 — 780 nm
  - CIE 1931 x,y and UCS 1976 u'v'
  - Tolerance ellipse and alarm
  - CCT and IRC
  - Dominant and peak wavelengths
- Data are exported or saved in Excel or text formats. The user manages to develop his own software application using the DLL controls library. Development is possible under various environment such as C++, Excel and LabView.



### Distributed in the UK and Ireland by


Photonics & Analytical Marketing Ltd  
7 Whitechapel Close  
Leeds, LS8 2PT  
Tel 0113 293 1985 Fax 0113 228 0204  
E-mail [sales@pandamarket.com](mailto:sales@pandamarket.com)  
Web site [www.pandamarket.com](http://www.pandamarket.com)



### Exchangeable measurement heads

Probe4Light is the world first PC-based spectrophotometer to feature both a miniature design and the performances of a lab class device.

For a maximum flexibility, Probe4Light may be configured with various measurement heads. The user may switch easily between them thanks to the screw-type binding. Other measurement heads may be manufactured upon the customer requests.



### Integrating sphere option

The integrating sphere has a 30 mm inner diameter and is coated with BaSO<sub>4</sub>. The sphere design enables to trap all the emitted light source and delivers a good light homogeneity.

The input port is 7 mm diameter and hence may accommodate most of the LED chips (Lumileds, Cree, Nichia, Osram, Seoul Semiconductors...)



The sphere should be in contact with the light source to collect all the emitted light.

Distributed in the UK and Ireland by

Photonics & Analytical Marketing Ltd  
7 Whitechapel Close  
Leeds, LS8 2PT  
Tel 0113 293 1985 Fax 0113 228 0204  
E-mail [sales@pandamarket.com](mailto:sales@pandamarket.com)  
Web site [www.pandamarket.com](http://www.pandamarket.com)



# Probe4Light

## Spectral Characterization of Light Sources

### Diffuser option

The diffuser is an opal material with a 20 mm diameter. The large photosensitive area delivers a good signal homogeneity.

The diffuser configuration is recommended for working on the light sources with a large emissive area, e.g. the backlights, the LEDs with collimated lenses and the CFL tubes.



The diffuser enables to work at various distances (irradiance set-up)



### SMA-905 adapter option

The SMA-905 adapter is a female terminal adapter for connecting an optical fiber. The device may thus be connected to any optical accessories such as a large integrating sphere or a cosine corrector probe.

Majantys proposes various accessories to configure Probe4Light for your applications. Please contact us.



SMA-905 connector with an optical fiber and a cosine corrector (example)



Distributed in the UK and Ireland by

Photonics & Analytical Marketing Ltd  
7 Whitechapel Close  
Leeds, LS8 2PT  
Tel 0113 293 1985 Fax 0113 228 0204  
E-mail [sales@pandamarket.com](mailto:sales@pandamarket.com)  
Web site [www.pandamarket.com](http://www.pandamarket.com)



# Probe4Light

## Spectral Characterization of Light Sources

	Probe4Light + integrating sphere	Probe4Light + SMA905 adapter + accessories	Probe4Light + diffuser (opal)	Probe4Light + other (please contact us)
Spectrum	✓	✓	✓	✓
Color, CCT, IRC	✓	✓	✓	✓
Irradiance* (W/m <sup>2</sup> ) Illuminance* (lux)		✓	✓	
Power* (W) Luminous flux* (lm)	✓			
Radiance* (W/m <sup>2</sup> /sr) Luminance* (cd/m <sup>2</sup> )				✓


(\*) Relative measurements only, for side-by-side sample comparison.

### Configuration summary

Probe4Light is delivered with its accessories and its software in a rugged case. Various accessories are available. Please contact us if a requested item is not listed underneath.



Spectrophotometer  
body



PhotonProbe software  
and DLL function library



SMA-905 adapter (to  
optical fiber)



Rugged case with foam



Integrating sphere - 30  
mm diameter



Integrating sphere - 80  
mm diameter



Diffuser (opal) - 20 mm  
diameter



Diffuser (opal) - 6 mm  
diameter




Distributed in the UK and Ireland by

Photonics & Analytical Marketing Ltd  
7 Whitechapel Close  
Leeds, LS8 2PT  
Tel 0113 293 1985 Fax 0113 228 0204  
E-mail [sales@pandamarket.com](mailto:sales@pandamarket.com)  
Web site [www.pandamarket.com](http://www.pandamarket.com)

# Probe4Light

## Spectral Characterization of Light Sources

Mechanical Properties	
Dimensions	130 mm length 55 mm outer diameter
Weight	320 gr.
Material	Anodized, stiff aluminium
Optical Properties	
Focal Length	20 mm
Aperture	f/2.7
Diffraction Grating	Aberration corrected type IV concave holographic grating
Entrance Aperture	100µm width – 500µm height
Resolution	10 nm FWHM with 100 µm width slit 5 nm as option — NEW
Stray Light	<0.1% at 415nm with RG630
Dispersion	100 nm/mm
Fiber Optic Connector	SMA905 – 0.22 numerical aperture optical fiber
Sensor specifications	
Detector	CMOS linear sensor
Detector Range	200 — 1100 nm
Pixels Resolution	Selectable 128,256,512,1024
Dynamic Range	71 dB
Absolute QE at peak	60% at 675 nm
A/D Resolution	12-Bit conversion
Integration Time	1 ms to 1 sec. (adjustable)
Features	Non Destructive Read Capable High sensitivity



System Performances	
Geometry	Integrating sphere Diffuser (opal)
Spectral Range	380 — 780 nm
Integration Time	1 ms to 3 sec.
Chromaticity accuracy x,y (*1)	0.005 RMS
Chromaticity repeatability x,y (*1)	+/- 0.0005
Dominant Wavelength Accuracy (*1)	+/- 0.5 nm
Dominant Wavelength Repeatability (*1)	+/- 0.5 nm
Computer	
Operating System	Windows XP and Vista NI LabView 8.2 driver
Interface (standard)	USB 2.0, one port
Interface (option)	Asynchronous Serial RS232/ RS485 Modbus — Ethernet — SDIO Other (please contact us)
Software	Windows software graphical user interface .NET controls available
Electronics	
Power Consumption	Max. 100 mA @ 5V Self powered, no external power supply
Input/Ouput	3 hardware ports available
Trigger	Yes

Note 1: Under calibration conditions coming along with the material.



Distributed in the UK and Ireland by

Photonics & Analytical Marketing Ltd  
7 Whitechapel Close  
Leeds, LS8 2PT  
Tel 0113 293 1985 Fax 0113 228 0204  
E-mail sales@pandamarket.com  
Web site www.pandamarket.com